

NEW PRODUCT!





SNOW **7** HD **Plus**

New Technology, New Experience

Key Specifications:

- Screen: 7-inch touch screen; screen resolution: 1280*720p
- Camera: 13M pixels; auto focus/tap to focus
- Magnification:
 - 2.5X ~ 19X (non-text-to-speech function)
 - 1.7X ~ 8X (real-time viewing mode under text-to-speech function)
- Color modes: Full color & 10 high contrast color modes
- Text-to-speech (optional)
- Viewing modes: Close & Distance
- Battery: Approximately 3.5 hours of continuous use;
approximately 3.5 hours charging time



Distance view



Close view