

# *Distance Camera* *for Snow 12*



**A great companion for our portable video magnifier Snow 12, this optional distance camera is helpful in the classroom or in the workplace.**



**● Plug and Play**

Connect the distance camera to Snow 12 and automatically switch to distance view mode.

**● Clear Display**

With a viewing distance up to 35 meters (almost 115 feet), objects or text on the other side of the room can be viewed.

**● Portable**

The distance camera is lightweight (11oz/335g) and easy to carry around.

## Specifications:

- Camera: 5 megapixels
- Connection: USB Type-C
- Dimensions: 5.3"×3.9"×9.6"/13.5cm×10cm×24.5cm