

ZOOMAX[®]

PRODUCT CATALOG

Your Vision Product Expert



Handheld Video Magnifiers

Snow 7 HD Plus

■ The perfect mix of power and portability



- 7-inch HD touch screen
- OCR and speech (20 languages)
- Two viewing modes: Distance view, Near view
- 2.5-19x Zoom and 10 color combinations
- Storage for up to 1,000 images

M5 HD Plus

■ Our most affordable HD solution



- 5-inch HD touch screen
- Dual HD cameras for distance view and close view
- 2.3-16x Zoom and 10 color combinations
- Adjustable reading line and masks
- Storage for up to 60 images

Snow 4.3

■ Our most popular pocket-size magnifier



- 4.3-inch TFT screen
- Lightweight and easy to use
- 2-16x Zoom and 10 color combinations
- Auto focus camera

Hot Seller

Snow 12

■ The only 12-inch portable magnifier solution



- 12-inch full HD touch screen
 - OCR and speech (20 languages)
 - 2.7-19x Zoom and 10 color combinations
 - Storage for up to 1,000 images
 - Supports JPG, PDF, RTF, and TXT file formats
 - Foldable Stand* for a bigger writing space
- *Single Snow 12 without Foldable Stand is available.

Desktop Video Magnifiers

Aurora HD

■ 24" widescreen foldable desktop video magnifier



- Adjustable 24-inch HD widescreen display
- Up to 70x magnification
- 16 color combinations
- Adjustable reading line and masks
- Moveable X-Y table

Panda

■ Most affordable desktop magnifier



- Adjustable 19-inch widescreen display
- Up to 70x magnification
- 16 color combinations
- Moveable X-Y table
- Available in standard or high definition

Warranty and Technical Service

Our products are built to very high standards in our ISO-9001 certified production facility, and registered as Class I Medical Devices with the U.S. FDA.

All Zoomax products are covered by a comprehensive 2-year factory warranty, including all parts and labor.

If you need help or support, send email to support@zoomaxusa.com.



About Zoomax

Zoomax Technology Co., Limited was established in 2011 and is a world leader in the development and manufacture of solutions for people with low vision. With our wholly owned subsidiary Zoomax (USA) Inc. in the Boston, MA area, two service facilities in California and Georgia, and a network of authorized distributors around the country and around the world, we are here to serve your needs today and in the future.

 www.zoomaxusa.com

 sales@zoomaxusa.com

