


ZOOMAX[®]



SNOW 12

SIMPLE AND POWERFUL

Put **SNOW 12** on the reading material and enjoy reading time now!



Single Snow 12 without Foldable Stand is available.



Simple & intuitive

Six simple buttons make it easy and intuitive for everyone to enjoy.



Superior image quality

Our best-ever camera produces the sharpest images with great details and lively colors.



Slim and modern

Less than an inch thick, the modern design is stylish and comfortable to hold.

Writing is possible and more comfortable than ever with a **foldable stand!**



Comfortable writing space

Bigger writing space for you to use for writing, doodling, signing and creating!



OCR and text-to-speech

Quickly and easily convert printed text into speech with over 20 languages.



See virtually everything with ease

Hold a pill bottle, printed material, recipe or homework under and read it.

Technical Specifications

Display	12-inch full HD touch screen
Full HD camera	13 megapixels
Magnification	2.7x to 19x or 2.5x to 19x w/stand
Color modes	Full color & 10 high-contrast color modes
OCR and text to speech (optional)	20+ languages
File formats	JPG, PDF, RTF, and TXT
Foldable stand	Yes
Working height under camera	14 cm (5.5 in.)
Distance viewing	Yes
Panning	Touch screen panning & analog stick panning
Connect to TV/monitor/iPhone	Yes
Continuous use time	2.5 hours
Charging time	2 hours
Weight	960 g (2.12 lbs.) or 2,240 g (4.9 lbs.) w/stand
Dimensions (H/W/D)	29x22x2.4 cm (12x9x1 in.)