


LUNA 6

We design it as you need.



Wireless charging



Comfortable reading with handle



Lightweight and portable



● Comfortable and Portable

Luna 6 handheld video magnifier is a pocket-sized, 6-inch magnifier that you can carry around and use everywhere. Offering a comfortable handle and unique wireless charge feature, Luna 6 is a practical gadget for low vision people on the go.

● Technical Specifications

- Magnification: 2.5x - 19x
- Screen: 5.5" TFT screen
- Camera: 13MP
- Focus: auto focus and refocus manually
- Color modes: full color and 10 high contrast color modes
- Foldable writing stand: yes
- Built-in handle: yes
- Battery life: more than 4h continuous use
- Battery: rechargeable lithium-ion battery
- Charging: both wired and wireless charging
- Weight: 210 g/0.46 lbs
- Dimensions: 160 x 80 x 11 mm/6.3 x 3.1 x 0.4 in